

GLUTEAL IMPLANTS

Silicone gel-filled gluteal implants available in two shapes and surfaces, with planar or concave back

POLYTECH Health & Aesthetics GmbH
Alzheimer Str. 32
64807 Dieburg
Germany

Tel +49 (0) 6071 98630
Fax +49 (0) 6071 986330
E-Mail info@polytechhealth.com
www.polytechhealth.com



2.2.14 2019-09 – © 2019 POLYTECH Health & Aesthetics GmbH, Germany – M0202014

POLYTECH

Gluteal Implants

CE 0483

Product Description

Silicone gel-filled gluteal implants for compensation of soft tissue deficiencies are available with two different surfaces:

- POLYsmooth™ (smooth)
- MESMO® (microtextured)

The surgeon can choose among oval and round shapes of various sizes as well as between a concave or planar implant posterior. Thus he is in a position to achieve the optimal result for either augmentation or reconstruction of the gluteal region in patients of any sex.



Information on Materials

The shell of all gluteal implants consists of a strong, resilient and highly elastic silicone elastomer equipped with a barrier layer. It is filled with a highly cross-linked silicone gel of firm consistency. The surface of the implants is either smooth, POLYsmooth™, or micro-textured, MESMO®. The arithmetic surface roughness of MESMO® implants is approx. 25 µm +/-2mm.

Production of the implants is based on a quality management according to ISO 13485. All materials used are of implantable medical quality.

Indications

- Augmentation
- Contour correction
- Reconstruction, e.g. after trauma
- Implant exchange

Important Information

Gluteal implants are delivered individually packed and sterile.

Attention: Please note the detailed instructions for use included with each product.

Gluteal Implants – Concave posterior

Silicone-gel filled, individually packed, sterile

Order N°	Vol. (ml)	Dimensions (mm)			Posterior
		A	B	C	

Oval, POLYsmooth™

10300-115 ¹	115	80	120	27	concave
10300-215 ¹	215	105	142	30	concave
10300-290	290	114	152	33	concave
10300-360	360	125	163	34	concave
10300-450 ¹	450	132	174	38	concave



Oval, MESMO®

15300-115 ¹	115	80	120	27	concave
15300-215	215	105	142	30	concave
15300-290	290	114	152	33	concave
15300-360	360	125	163	34	concave
15300-450 ¹	450	132	174	38	concave



Round, POLYsmooth™

10301-200 ¹	200	107	107	35	concave
10301-250 ¹	250	112	112	41	concave
10301-300	300	118	118	44	concave
10301-360	360	127	127	43	concave
10301-500 ¹	500	141	141	49	concave



Gluteal Implants — Planar posterior

Silicone-gel filled, individually packed, sterile

Order N°	Vol. (ml)	Dimensions (mm)			Posterior
		A	B	C	

Oval, POLYsmooth™

10300-240 ¹	240	90	130	40	planar
10300-305	305	95	139	45	planar
10300-365	365	100	148	48	planar
10300-420	420	109	152	50	planar



Oval, MESMO®

15300-240	240	90	130	40	planar
15300-305	305	95	139	45	planar
15300-365	365	100	148	48	planar
15300-420	420	109	152	50	planar



Round, POLYsmooth™

10301-225 ¹	225	107	107	35	planar
10301-275	275	112	112	41	planar
10301-340	340	118	118	44	planar
10301-390	390	128	128	44	planar
10301-560	560	141	141	49	planar



Approximate measurements

¹ Supply after order ≈ 8-10 weeks