



Testicular Implants

Filled with silicone gel or soft silicone elastomer, with reinforced cap or needle guide, each type available in five sizes



With reinforced cap



With needle guide

Implants made by POLYTECH – Quality made in Germany



Testicular Implants, € 0483

Product Description

Testicular implants are available in two different types, both suitable for attachment to the tissue.

The implants of the type filled with silicone gel are fitted with a short piece of removable teflon to guide the needle and prevent a puncture of the implant.

The implants of the type filled with soft silicone elastomer have a reinforced cap that allows the attachment of the implant with a suture through the cap.

Both types are available in five different sizes.

Information on Materials

The multi-layer shell of the testicular implants is equipped with a diffusion barrier and has a smooth POLYsmooth™ surface. The implants are filled with either a highly cross-linked silicone gel or a soft silicone elastomer.

Production of the implants is based on a quality management according to ISO 9001 and ISO 13485. All materials used are of implantable medical quality.

Indications

- Reconstruction after orchidectomy
- Reconstruction after trauma
- Correction of asymmetry
- Congenital defects
- Implant replacement


Testicular implants

Order n°	Volume (ml)	Dimensions (cm)	
		A	B

with needle guide, POLYsmooth™

Filled with silicone-gel, smooth surface, individually packed, sterile, € 0483


10100-005	5	2.0	2.2
10100-012	12	2.6	3.3
10100-021	21	3.2	4.2
10100-027	27	3.3	4.8
10100-036	36	3.9	5.0



with reinforced cap, POLYsmooth™

Filled with soft silicone elastomer, smooth surface, individually packed, sterile, € 0483

10101-005	5	2.0	2.2
10101-012	12	2.6	3.3
10101-021	21	3.2	4.2
10101-027 ¹	27	3.3	4.8
10101-036 ¹	36	3.9	5.0



approximate measurements

On stock	¹ Supply after order ≈ 6 weeks
----------	---

Important Information

Testicular implants are delivered individually packed and sterile.

CAUTION: The teflon needle guide of the silicone gel-filled implants must always be removed before wound closure!

Attention: Please note the detailed instructions for use included with each product.



POLYTECH Health & Aesthetics GmbH

Altheimer Str. 32 • 64807 Dieburg • Germany

+49 (0) 60 71 98 63 0 • +49 (0) 60 71 98 63 30

info@polytechhealth.com

www.polytech-health-aesthetics.com