



Bustier without cups white  
(with cups see left)



Bustier with cups black



Bustier without cups black



Bustier combined with bandage

## **POLYTECH Bustiers and Compressa™ Bandages**

post-operative compression bras  
and bandages for optimising  
breast-surgery results

Products made by POLYTECH – QUALITY made in Germany



# POLYTECH Bustiers and Compressa™ Bandages, €€

## Product Description



Bustier back view

POLYTECH Bustiers are available in seven sizes, with or without cups, and in the colours black and white. They provide an optimal postoperative care after mammary augmentation, reconstruction or reduction.

The bustiers are made of mixed fibres with a high cotton content for greater patient comfort. They can be washed by hand at 30°C. Adjustable velcros, a zip-front opening and three hook-and-eye fasteners guarantee an easy handling.

The Compressa™ Breast Bandage for post-operative compression therapy is available in two sizes: 90cm length (medium) and 95cm length (large), both with a width of 6.5cm. The bandage is easy to apply: its size and pressure are adapted flexibly via the Velcro® fastener that can be fastened anywhere along the bandage.

The bandage material – 89% polyamide, 11% elastane – is highly elastic, soft and well tolerable to the skin, thus very comfortable. The bandages may be machine-washed at 40°C.

## Advantages

- Reduced risk of the occurrence of post-operative hematoma, oedema, and hypertrophic scars
- Fixation of implants and wound-care bandages in their position
- Facilitation of wound inspection and handling

## Important Information

POLYTECH Bustiers and Compressa™ Bandages are provided individually packed and non-sterile.

**Attention:** Please note the detailed instructions for use included with each product.

## POLYTECH Bustiers

Post-operative compression bras, individually packed, non sterile, €€

Order n°	Under-bust measurement (cm)	Size
<b>Bustier with cups, WHITE</b>		
2110 P	< 84	Petite
2110 S	84 – 89	Small
2110 M	89 – 94	Medium
2110 L	94 – 104	Large
2110 XL	104 – 114	eXtra Large
2110 XXL	114 – 124	XX Large
2110 XXXL	124 – 134	XXX Large



WHITE

Order n°	Under-bust measurement (cm)	Size
<b>Bustier with cups, BLACK</b>		
2120 P	< 84	Petite
2120 S	84 – 89	Small
2120 M	89 – 94	Medium
2120 L	94 – 104	Large
2120 XL	104 – 114	eXtra Large
2120 XXL	114 – 124	XX Large
2120 XXXL	124 – 134	XXX Large



BLACK

Order n°	Under-bust measurement (cm)	Size
<b>Bustier without cups, WHITE</b>		
2111 P	< 84	Petite
2111 S	84 – 89	Small
2111 M	89 – 94	Medium
2111 L	94 – 104	Large
2111 XL	104 – 114	eXtra Large
2111 XXL	114 – 124	XX Large
2111 XXXL	124 – 134	XXX Large



WHITE

Order n°	Under-bust measurement (cm)	Size
<b>Bustier without cups, BLACK</b>		
2121 P	< 84	Petite
2121 S	84 – 89	Small
2121 M	89 – 94	Medium
2121 L	94 – 104	Large
2121 XL	104 – 114	eXtra Large
2121 XXL	114 – 124	XX Large
2121 XXXL	124 – 134	XXX Large



BLACK

## Compressa™ Bandage

Post-operative compression bandage, individually packed, non sterile, €€

Order n°	Specification
2221M	MEDIUM, black, length 90cm, €€
2221L	LARGE, black, length 95cm, €€



Bandage side view & Velcro® fastening



POLYTECH Health & Aesthetics GmbH

Altheimer Str. 32 • 64807 Dieburg • Germany

☎ +49 (0) 60 71 98 63 0 • 📠 +49 (0) 60 71 98 63 30

✉ info@polytechhealth.com • www.polytechhealth.com